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Read the instructions for use carefully before starting to use the product.
Keep these instructions, as you may have to read them again.

Sibokit

Breath test for the diagnosis of bacterial overgrowth

1. Manufacturer's information

Registered corporate name and address

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2. Case contents

The SIBO-KIT package contains:

- 1 sachet of lactitol monohydrate (10 g)
- 12 breath sample collection tubes.
- 2 flexible tubes.
- 1 instruction for use leaflet.
- Sample identification labels.

3. Product type

SIBO-KIT is an in vitro diagnostic test for bacterial overgrowth. The product is used for taking breath samples from patients who have symptoms associated with bacterial overgrowth.

4. SIBO-KIT conservation

Store in a dry and cool place.

5. Special instructions for use

It is important that these instructions be followed precisely to ensure the reliability and reproducibility of the diagnosis. It is recommended that the test be performed in the presence of qualified personnel.

PREPARING THE PATIENT BEFORE THE TEST

- The conduct of this test is not recommended within 4 weeks of treatment with antibiotics, osmotic laxatives, and/or cleansing enemas, or for patients who have undergone endoscopic or radiological testing requiring bowel cleansing.
- It is advised that the consumption of probiotics be discontinued 2 weeks before the test.
- A diet low in fermentable carbohydrates should be consumed for 24 hours prior to the test to achieve optimal baseline levels of H₂ and/or CH₄.
- The patient must fast for a minimum period of 12 hours.
- On the day of the test, the patient shall wash their mouth exclusively by rinsing it with a sorbitol-free chlorhexidine solution.
- The patient should refrain from chewing gum, smoking, or performing any physical effort starting at least 30 minutes before the start of the test and until its completion.
- Only non-carbonated water can be consumed during the test.

TEST PROCEDURE

1. The breath collection tubes must be labeled with the sampling times as shown in the following table:

Tube no.	1	2	3	4	5	6
Time	BASELINE (0 min)	POST (15 min)	POST (30 min)	POST (45 min)	POST (60 min)	POST (75 min)

Tube no.	7	8	9	10	11	12
Time	POST (90 min)	POST (105 min)	POST (120 min)	POST (135 min)	POST (150 min)	POST (165 min)

2. Before administering the substrate, the patient must **blow through a flexible tube** located at the bottom of **tube 1** in order to obtain the **BASELINE** breath sample.
*It is recommended that the patient hold their breath for approximately 30 seconds so that the concentration of H₂, CH₄, and CO₂ in their breath is maximized. The sample must contain the last portion of the exhaled air. The glass tube should be capped quickly while the patient continues to blow in order to avoid sample losses.
3. The lactitol substrate should be then administered. For adult patients, the recommended dose is 0.5 g/kg of weight (up to a maximum dose of 10 g) dissolved in a glass containing 250 ml of water. For pediatric patients, the substrate must be dissolved in 125 ml of water. The test can also be performed using glucose at a dose of 1g/kg of weight (up to a maximum dose of 100 g) or lactulose at a dose of 0.5 g/kg of weight (up to a maximum dose of 10 g) as substrates, both dissolved in a glass with 250 ml water. Identify the sample by sticking the label on the front of the box included in the Kit corresponding to the substrate used.
4. After administering the substrate, breath samples should be collected into the appropriate **POST tubes** every 15 minutes according to the table above. The sampling procedure shall be the same as described in point 2.
5. Identify the kit by sticking the label with symbol ① to the corresponding area of the case. The label with the word "**PATIENT**" must be kept by the patient, and the one identified with the word "**CUSTOMER**" must be kept by the site as proof.
6. Send the kit together with the **PATIENT'S SYMPTOMS TABLE** to the Analysis Service of **ISOMED PHARMA S.L.**

TEST PROTOCOL

It is recommended that the test be performed in the presence of qualified personnel

SIBO-KIT Protocol

1

Collect the **pre-dosing** sample (BASELINE).

2

Take the substrate dissolved in 250/125 ml of water.

3

Wait **15 minutes**.

4

Collect a **post-dosing** sample every **15 min**, for a total of 165 min.

5

Identification of samples and **shipping** for analysis.

6. Incompatibilities

Not applicable due to the product's characteristics.

7. Warnings

This product contains **lactitol**. Patients with hereditary fructose or galactose intolerance, galactosemia, glucose-galactose malabsorption, intestinal obstruction, or serum electrolyte imbalance should not use this product.

If for any reason you need to repeat the test, please note that at least 24 hours must elapse between the first and second administration of the lactitol substrate.

This medical device is **for single use only**. Cross-contamination may occur if the product is used for more than one patient, thus rendering the test useless.

The sample collection tubes must be closed properly without exerting excessive pressure. It is recommended that the cap be screwed on until it no longer turns easily, and to then carefully turn it another quarter turn.

8. Reporting of incidents

If the patient experiences any type of adverse incident, please contact a qualified professional. You may also report any incidents directly to the manufacturer through the following email: farmacovigilancia@isomed.com

INSTRUCTIONS FOR SENDING BREATH SAMPLES FOR ANALYSIS

1.- Fill in the following data from the label identified by symbol "Ⓢ" that you will find inside the case:

- Source code
- Age
- Sex (please mark with an X)
- Amount of substrate

2.- Once the label has been filled in, stick it to the area of the case with symbol "Ⓢ" reserved for this purpose.

3.- Save the rest of the labels identifying the sample ("CUSTOMER/PATIENT") as proof until receipt of the result.

4.- Send the kit (case containing the tubes after the breath collection, instructions for use with the patient's symptoms table filled in, and label identified by symbol "Ⓢ" attached to the corresponding area of the case) to the following address:

Isomed Pharma, S.L.
Avenida Alberto Alcocer, nº 46A, 4ºB
28016 Madrid
Email: analisis@isomed.com

It is very important that the following sender's data be indicated in each sample for the purpose of subsequently resending the result.

- **Address**
- **Email address (if any)**
- **Contact telephone number**
- **List of samples sent**

The following form was prepared to help health professionals assess the result of the Breath Test and diagnose bacterial overgrowth based on a patient's symptoms.

IMPORTANT NOTE: Once the form has been completed by the patient, please insert it into the kit and send it to the Analysis Service of ISOMED PHARMA S.L.

PATIENT'S INFORMATION

Source code: _____

Substrate administered: _____

Amount of substrate: _____

Age: _____

Sex: _____

Male

Female

ASSESSMENT OF THE PATIENT'S SYMPTOMS

Mark the symptoms that the patient experienced during the test with an "X".

SYMPTOMS	BEFORE (Previous fast)	DURING (Sampling)	AFTER (Test end)
Abdominal pain or cramping			
Gases (belching or flatulence)			
Bowel or abdominal sounds			
Nausea or vomiting			
Constipation			
Diarrhea			
Abdominal inflammation or swelling			
Heartburn			
Dyspepsia or digestive disorders			
Other symptoms (please specify):			

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Lactokit

Breath test for the diagnosis of lactose intolerance

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2. Case contents

The LACTO-KIT package contains:

- 1 sachet of lactose monohydrate (25g)
- 8 breath sample collection tubes
- 2 flexible tubes
- 1 instruction for use leaflet
- Sample identification labels

3. Product type

LACTO-KIT is an in vitro diagnostic test for lactose malabsorption or intolerance. The product is used for taking breath samples from patients who have symptoms associated with lactose intolerance.

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2. Before administering the substrate, the patient must blow through a **flexible tube** located at the bottom of **tube1** in order to obtain the **BASELINE** breath sample.
*It is recommended that the patient hold their breath for approximately 30 seconds so that the concentration of H₂, CH₄, and CO₂ in their breath is maximized. The sample must contain the last portion of the exhaled air. The glass tube should be capped quickly while the patient continues to blow in order to avoid sample losses.
3. The **lactose** substrate should then be administered. For adult patients, the recommended dose is 1 g/kg of weight (up to a maximum dose of 25 g) dissolved in a glass containing 250 ml of water. For pediatric patients, the substrate must be dissolved in 125 ml of water.
*See section 7 (Warnings) in case of lactose intolerance.
4. After administering the substrate, breath samples should be collected into the appropriate **POST** tubes every 25 minutes according to the table above. The sampling procedure shall be the same as described in point 2.
5. Identify the kit by sticking the label with symbol ① to the corresponding area of the case. The label with the word "**PATIENT**" must be kept by the patient, and the one identified with the word "**CUSTOMER**" must be kept by the site as proof.
6. Send the kit together with the **PATIENT'S SYMPTOMS TABLE** to the Analysis Service of **ISOMED PHARMA S.L.**

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LACTO-KIT Protocol

1 Collect the **pre-dosing** sample (BASELINE).

2 Take the substrate dissolved in 250/125 ml of water.

3 Wait **25 minutes**.

4 Collect a **post-dosing** sample every **25 min**, for a total of 175 min.

5 **Identification** of samples and **shipping** for analysis.

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7. Warnings

This product contains **Lactose**. Patients with hereditary fructose or galactose intolerance, galactosemia, glucose-galactose malabsorption, intestinal obstruction, or serum electrolyte imbalance should not use this product.

May contain traces of milk. Patients with allergy to milk protein should not use this product.

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8. Reporting of incidents

If the patient experiences any type of adverse incident, please contact a qualified professional. You may also report any incidents directly to the manufacturer through the following email: **farmacovigilancia@isomed.com**

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1.- Fill in the following data from the label identified by symbol "⊕" that you will find inside the case:

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2.- Once the label has been filled in, stick it to the area of the case with symbol "⊕" reserved for this purpose.

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4.- Send the kit (case containing the tubes after the breath collection, instructions for use with the patient's symptoms table filled in, and label identified by symbol "⊕" attached to the corresponding area of the case) to the following address:

Isomed Pharma, S.L.
Avenida Alberto Alcocer, nº 46A, 4ºB
28016 Madrid, Spain
Email: analisis@isomed.com

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The following form was prepared to help health professionals assess the result of the Breath Test and diagnose a lactose intolerance based on a patient's symptoms.

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TIME SYMPTOMS	BEFORE (fasted)	DURING 25min	DURING 50min	DURING 75min	DURING 100min	DURING 125min	DURING 150min	DURING 175min	AFTER the test
	Abdominal pain or cramping								
Gases (belching or flatulence)									
Bowel or abdominal sounds									
Nausea or vomiting									
Constipation									
Diarrhea									
Abdominal inflammation or swelling									
Heartburn									
Dyspepsia or digestive disorders									
Other symptoms (please specify):									

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2. Case contents

The FRUCTO-KIT package contains:

- 1 sachet of granulated crystalline fructose monohydrate (25g)
- 8 breath sample collection tubes
- 2 flexible tubes
- 1 instruction for use leaflet
- Sample identification labels

3. Product type

FRUCTO-KIT is an in vitro diagnostic test for fructose malabsorption or intolerance. The product is used for taking breath samples from patients who have symptoms associated with fructose intolerance.

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2. Before administering the substrate, the patient must blow through a flexible tube located at the bottom of tube1 in order to obtain the **BASELINE** breath sample.
*It is recommended that the patient hold their breath for approximately 30 seconds so that the concentration of H₂, CH₄, and CO₂ in their breath is maximized. The sample must contain the last portion of the exhaled air. The glass tube should be capped quickly while the patient continues to blow in order to avoid sample losses.
3. The fructose substrate should then be administered. For adult patients, the recommended dose is 1 g/kg of weight (up to a maximum dose of 25 g) dissolved in a glass containing 250 ml of water. For pediatric patients, the substrate must be dissolved in 125 ml of water.
In case intolerance to fructose + sorbitol is to be evaluated, dissolve 0,50 g sorbitol/kg of weight, up to a maximum of 5 g in the same glass and identify the sample by sticking the label "FRUCTOSE+SORBITOL" (included in the Kit) on the front of the box.
*See section 7 (Warnings) in case of fructose intolerance.
4. After administering the substrate, breath samples should be collected into the appropriate **POST tubes** every 25 minutes according to the table above. The sampling procedure shall be the same as described in point 2.
5. Identify the kit by sticking the label with symbol ① to the corresponding area of the case. The label with the word "**PATIENT**" must be kept by the patient, and the one identified with the word "**CUSTOMER**" must be kept by the site as proof.
6. Send the kit together with the **PATIENT'S SYMPTOMS TABLE** to the Analysis Service of **ISOMED PHARMA S.L.**

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Wait **25 minutes**.

4

Collect a **post-dosing** sample every **25 min**, for a total of 175 min.

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The following form was prepared to help health professionals assess the result of the Breath Test and diagnose a fructose intolerance based on a patient's symptoms.

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TIME	BEFORE (tasted)	DURING 25min	DURING 50min	DURING 75min	DURING 100min	DURING 125min	DURING 150min	DURING 175min	AFTER the test
SYMPTOMS									
Abdominal pain or cramping									
Gases (belching or flatulence)									
Bowel or abdominal sounds									
Nausea or vomiting									
Constipation									
Diarrhea									
Abdominal inflammation or swelling									
Heartburn									
Dyspepsia or digestive disorders									
Other symptoms (please specify):									

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7. Warnings

This product contains **fructose**. Patients with hereditary fructose intolerance should not take this product.

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