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Read carefully the instructions for use before starting to use the product.
Keep these instructions, as you may have to read them again.

DisbioKit

Intestine

Kit for stool sample collection test for the diagnosis of intestinal dysbiosis.

1. Manufacturer's information

Registered corporate name and address

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www.isomed.com

2. Case contents

The DISBIO-KIT INTESTINE package contains:

- 1 Sample collection paper.
- 1 sample container.
- 1 Glove.
- 1 UN3373 bag.
- 1 instruction leaflet.
- Identification labels.

3. Product type

DISBIO-KIT INTESTINE is a kit for stool sample collection for the diagnosis of intestinal dysbiosis. The product is intended for the collection of stool samples in order to analyze the intestinal microbiota, as well as other other specific parameters related to intestinal health and several disorders.

4. Conservation of DISBIO-KIT INTESTINE

Store in a cool and dry place until use.

Once the sample is collected, store at a temperature of 2-8°C.

5. Special instructions for use



It is important to follow these instructions precisely to ensure the reliability and reproducibility of the diagnosis.

PREPARATION FOR THE TEST

No exceptional preparation is required. Maintain your usual dietary habits and collect the stool sample anytime during the day.

It is not necessary to take any previous medication or nutritional supplement. Please record any medications and nutritional supplements you are currently taking, as some may change the characteristics of faeces.

TEST PROCEDURE

 <p>1</p>	<p>Clean and dry the surface of the toilet.</p>	 <p>6</p>	<p>Collect a stool sample with the spatula in the container. Repeat this operation 3 times at different points of the sample.</p>
 <p>2</p>	<p>Unfold the sample collection paper and stick it on the edge of the toilet.</p>	 <p>7</p>	<p>Introduce the sample into the collection container with the single amount filled in the spatula of the container.</p>
 <p>3</p>	<p>Make sure the sample collection paper has the shape of a bowl.</p>	 <p>8</p>	<p>Screw the lid to seal the containers and stick the label with the text "RECIPIENT".</p>
 <p>4</p>	<p>Proceed to make your bowel movements.</p>	 <p>9</p>	<p>Put the container in the UN3373 plastic bag. Introduce the sample in the freezer at 2-8°C until it is picked up.</p>
 <p>5</p>	<p>Wear the disposable glove.</p>	 <p>x 2-3</p> <p>10</p>	<p>Dispose the sample collection paper in the toilet and flush.</p>

Identify the kit sticking the label with the text CAJA/BOX in the corresponding area of the case. The label with the word "PATIENT" must be preserved by the patient, and the label identified with the word "VOLANTE" must be adhered to the analysis request flyer.

Send the kit together with the ANALYSIS REQUEST FLYER to ISOMED DIAGNOSTICS Analysis Laboratory.

6. Incompatibilities

Not applicable due to the product's characteristics.

7. Warnings

This medical device is for single use only. Cross-contamination may occur if the product is used by more than one patient, thus rendering the test useless.

Dispose of all samples and materials used as biohazard waste

8. Reporting of incidents

If the patient experiences any type of adverse incident, please contact a qualified professional. You may also report any incidents directly to the manufacturer through the following email: farmacovigilancia@isomed.com

INSTRUCTIONS FOR SENDING STOOL SAMPLES FOR ANALYSIS

- 1.- Stick the label "VOLANTE" on the analysis request form.
- 2.- Stick the CAJA/BOX label to the case.
- 3.- Keep the "PACIENTE" label as proof until the reception of the results report.
- 4.- Put the UN3373 plastic bag with the sample container inside the case, along with the filled analysis request form. Place the label identified by the text CAJA/BOX in the corresponding area of the case.
- 5.- To send the sample, please contact ISOMED DIAGNOSTICS S.L. by calling: 911354048 to arrange the collection of the kit.