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Read the instructions for use carefully before starting to use the product.
Keep these instructions, as you may have to read them again.

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Kit for breath sample collection for the diagnosis of sugars intolerance

1. Manufacturer's information

Registered corporate name and address

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- 8 breath sample collection tubes
- 2 flexible tubes
- 1 instruction for use leaflet
- Sample identification labels

3. Product type

PRUEBAS FUNCIONALES is a kit for breath sample collection for the diagnosis of sugars intolerance or malabsorption. The product is used for taking breath samples from patients who have symptoms associated with sugar intolerance or malabsorption.

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No special storage conditions required.

5. Special instructions for use

It is important that these instructions be followed precisely to ensure the reliability and reproducibility of the diagnosis. It is recommended that the test be performed in the presence of qualified personnel.

PREPARING THE PATIENT BEFORE THE TEST

- The conduct of this test is not recommended within 4 weeks of treatment with antibiotics, osmotic laxatives, and/or cleansing enemas, or for patients who have undergone endoscopic or radiological testing requiring bowel cleansing.
- It is advised that the consumption of probiotics be discontinued 2 weeks before the test.
- A diet low in fermentable carbohydrates should be consumed for 24 hours prior to the test to achieve optimal baseline levels of H₂ and/or CH₄.
- The patient must fast for a minimum period of 8 hours.
- On the day of the test, the patient shall wash their mouth exclusively by rinsing it with a sorbitol-free chlorhexidine solution.
- The patient should refrain from chewing gum, smoking, or performing any physical effort starting at least 30 minutes before the start of the test and until its completion.
- Only non-carbonated water can be consumed during the test.

TEST PROCEDURE

1. The breath collection tubes must be labeled with the sampling times as shown in the following table:

Tube no	Tube 1	Tube 2	Tube 3	Tube 4	Tube 5	Tube 6	Tube 7	Tube 8
Time	BASELINE (0 min)	POST (25 min)	POST (50 min)	POST (75 min)	POST (100 min)	POST (125 min)	POST (150 min)	POST (175 min)

2. Before administering the substrate, the patient must **blow through a flexible tube** located at the bottom of **tube1** in order to obtain the **BASELINE** breath sample.

*It is recommended that the patient hold their breath for approximately 30 seconds so that the concentration of H₂, CH₄, and CO₂ in their breath is maximized. The sample must contain the last portion of the exhaled air. The glass tube should be capped quickly while the patient continues to blow in order to avoid sample losses.

3. The substrate of choice is then administered dissolved in a glass with 250 ml of water if the patient is an adult. For pediatric patients, the substrate must be dissolved in 125 ml of water. The recommended dose of each substrate is:

D-xylose	2 g/kg of weight up to a maximum of 25 g.
Sucrose	2 g/kg of weight up to a maximum of 50 g.
Galactose	1 g/kg of weight up to a maximum of 25 g.
Maltose	1 g/kg of weight up to a maximum of 25 g.
Sorbitol	0,5 g/kg of weight up to a maximum of 10 g.

4. After administering the substrate, breath samples should be collected into the appropriate **7 POST tubes** every 25 minutes according to the table above. The sampling procedure shall be the same as described in point 2.

5. Identify the kit by sticking the label with the text **CAJA/BOX** to the corresponding area of the case. The label with the word **"PATIENT"** must be kept by the patient, and the one identified with the word **"CUSTOMER"** must be kept by the site as proof.

Indicate the **substrate and dose** administered for the test in the corresponding box of the package.

6. Send the kit together with the **PATIENT'S SYMPTOMS TABLE** to the Analysis Service of **ISOMED PHARMA S.L.**

TEST PROTOCOL

It is recommended that the test be performed in the presence of qualified personnel

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1 Collect the **pre-dosing** sample (BASELINE).

2 Take the substrate dissolved in 250/125 ml of water.

3 Wait **25 minutes**.

4 Collect a **post-dosing** sample into each of the **7 tubes every 25 min**,

5 **Identification** of samples and **shipping** for analysis.

7. Incompatibilities

Not applicable due to the product's characteristics.

8. Warnings

This medical device is **for single use only**. Cross-contamination may occur if the product is used for more than one patient, thus rendering the test useless.

The sample collection tubes must be closed properly without exerting excessive pressure. It is recommended that the cap be screwed on until it no longer turns easily, and to then carefully turn it another quarter turn.

Dispose of all samples and materials used as biohazard waste.

9. Reporting of incidents

If the patient experiences any type of adverse incident, please contact a qualified professional. You may also report any incidents directly to the manufacturer through the following email:

farmacovigilancia@isomed.com.

INSTRUCTIONS FOR SENDING BREATH SAMPLES FOR ANALYSIS

- 1.- Fill in the following data from the label identified by the text **CAJA/BOX** that you will find inside the case:
 - Source code
 - Age
 - Sex (please mark with an X)
 - Administered substrate
 - Amount of substrate
- 2.- Once the label has been filled in, stick it to the area of the case with the text **CAJA/BOX** reserved for this purpose.
- 3.- Save the rest of the labels identifying the sample ("CUSTOMER/PATIENT") as proof until receipt of the result.
- 4.- Send the kit (case containing the tubes after the breath collection, instructions for use with the patient's symptoms table filled in, and label identified by the text **CAJA/BOX** attached to the corresponding area of the case) to the following address:

Isomed Diagnostics
Calle de Julian Camarillo, 7, 3^oB
28037 Madrid
Email: analisis@isomed.com

It is very important that the following sender's data be indicated in each sample for the purpose of subsequently resending the result.

- **Address**
- **Email address (if any)**
- **Contact telephone number**
- **List of samples sent**

The following form was prepared to help health professionals assess the result of the Breath Test and diagnose a sugar intolerance based on a patient's symptoms.

IMPORTANT NOTE: Once the form has been completed by the patient, please insert it into the kit and send it to the Analysis Service of ISOMED PHARMA S.L.

PATIENT'S INFORMATION

Source code: _____

Substrate administered: _____

Amount of substrate _____

Age: _____

Sex: Male Female

ASSESSMENT OF THE PATIENT'S SYMPTOMS

Mark the symptoms that the patient experienced during the test with an "X".

TIEMPO SÍNTOMAS	BEFORE (fasted)	DURING 25min	DURING 50min	DURING 75min	DURING 100min	DURING 125min	DURING 150min	DURING 175min	AFTER the test
Abdominal pain or cramping									
Gases (belching or flatulence)									
Bowel or abdominal sounds									
Nausea or vomiting									
Constipation									
Diarrhea									
Abdominal inflammation or swelling									
Heartburn									
Dyspepsia or digestive disorders									
Other symptoms (please specify):									

Comments: _____